

### Luxir

#### **Pure MSC Exosomes in normal saline with high-potency Amino Acids**

Luxir is a pure MSC exosome product for topical use. The product contains proteins, enzymes, growth factors, and exosomes derived from mesenchymal stem cells (MSCs). Luxir is an isolated, cell-free, placental mesenchymal stem cell-derived exosome re-suspended in 0.9% normal saline; every lot is tested for sterility and endotoxin and is quantified AND qualified via protein concentration, dSTORM visualization, NanoSight sizing, mRNA sequencing, Mass Spectrometric protein sequencing, flow cytometric surface marker testing, and visual purity assay for the world's first reproducible pure exosome product manufactured to FDA drug and vaccine production guidelines.

Name	Volume	Ref	Biologically Active Exosomes	Price 1 - 5	Price 6 - 14	Price 15+
Luxir	1.0 mL	LXR0101	1 Trillion Microvesicles	\$350	\$300	\$250
Luxir	5.0mL	LXR0501	5 Trillion Microvesicles	\$1,200	\$1,100	\$1,000
Luxir+	1.0mL	LXRP0101+	3 Trillion Microvesicles	\$800	\$700	\$600
Luxir+	5.0mL	LXRP0501+	15 Trillion Microvesicles	\$3,000	\$2,750	\$2,500

Exosome numbers stated within are for Kimera potency purposes **only**.

**Do not compare exosome numbers of Kimera products to other products derived from alternate cell sources.**

Not a tissue product as defined by FDA guidelines.

### Vive

#### **MSC Exosome Growth Factor Skincare Serum with Hyaluronic Acid, Vitamin C, High-potency Amino Acids, Rose Water, plus Trehalose.**

Vive is a topical solution specially developed for skincare purposes. The product contains proteins, enzymes, growth factors, and exosomes derived from mesenchymal stem cells (MSCs). Vive is an isolated, cell-free, placental mesenchymal stem cell-derived exosome re-suspended in 0.9% normal saline.

Name	Volume	Ref	Biologically Active Exosomes	Price	Price 6 - 14	Price 15+
Vive	5.0 mL	VPI0501	2 Trillion Microvesicles	\$135	\$135	\$135

### **SHIPMENT AND STORAGE OF KIMERA PRODUCTS**

- Kimera Labs' exosomes stored at -80°C have a two-year or one-year shelf life if stored at -20C.
- **DO NOT STORE ANY KIMERA PRODUCT IN LIQUID NITROGEN OR A CRYOGENIC FREEZER.** Kimera products are cell-free; this is unnecessary. Kimera Labs may not be responsible if vials fracture with improper storage.
- All Kimera Labs products are shipped FedEx Priority Plus Overnight on dry ice, scheduled to arrive by noon (unless First AM Delivery is requested). Tracking information will be provided.
- Be prepared to receive – a *signature is required for delivery*. Place the product directly in the freezer.
- Keep product frozen until use; thaw at room temperature. Use thawed product within 24 hours (kept refrigerated), do not refreeze.

**Send the completed order form to [ORDERS@KIMERALABS.COM](mailto:ORDERS@KIMERALABS.COM) or**

**Order online at <https://store.kimeralabs.net>**